

## ROSETT GLASS COOLING CELLS

The Rosett cooling cell enables uniform treatment at low temperatures. The cell is placed in an ice bath. The ultrasonic energy forces the sample to circulate repeatedly under the probe and through-out the cooling arms.

30 ml Rosett cooling cell.  
Part No. 830-00003



## GLASS COOLING CELLS\*

10 ml cooling cell with water jacket.  
Part No. 830-00009

100 ml cooling cell with water jacket.  
Part No. 830-00010



\*Inlet and outlet require  $\frac{3}{8}$ " (9.5 mm) I.D. tubing.